

## RESISTANCE OF POLYETHYLENE TO CHEMICAL PRODUCTS

### FEATURES OF POLYETHYLENE PE / POLYPROPYLENE PP

- High resistance to chemical products and sunlight.
- High resistance to impacts.
- 100% recyclable.

### RESISTANCE COMPARATIVE

	PE / PP	ABS
ALCOHOL	High resistance	Medium resistance
CHLORHEXIDINE	High resistance	Medium resistance
IODINE (BETADINE)	High resistance	Medium resistance
ACETONE	Medium resistance	Does not resist



PE WORK SURFACE



PE FOLD DOWN TRAY



TROLLEY 480 x 480 mm  
PP DRAWERS, 3 COLOURS



TROLLEY 640 x 480 mm / 900 x 630 mm  
PE ISO DRAWER FRONTS, 6 COLOURS

INSAUSTI manufactures the drawer fronts in polyethylene since 1997 and the work surface since 2013. In this period no claims have been received for damages caused by chemical products.