

## TECHNICAL FEATURES

### SQM SERIES SQUARE TROLLEY

480 x 480 mm

#### EXTENSIVE MOBILITY TROLLEYS

#### WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>



#### SQUARE DRAWER

- Polypropylene drawers, PP.
- Resistance to chemicals similar to the polyethylene.
- Manufactured in a single piece with slots for optional dividers.
- Heights 57, 130 mm.
- Colour white, blue, red.
- Total opening, steel ball-bearing slides with retention system.
- Front stop, it protects the slides against strong impacts when closing.
- Identification cards ISO norm 7810.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>



Elevated push handle

Trolley prepared to mount on accessories at destination

Monobloc structure of AISI-304 stainless steel

Polyamide handle, PA, ergonomic  
Easy-to-clean

**SIMPLE WHEELS**  
Ø 125 mm

2 with brake

2 Antistatic

