

## TECHNICAL FEATURES

### LOW SIMPLE TROLLEY, REDUCED HEIGHT

640 x 480 mm

#### STATIC TROLLEYS

#### WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>



#### ISO DRAWERS 600 x 400 mm

- ABS ISO drawers of standard dimensions according to ISO 3394, compatible and interchangeable with ISO furniture of hospitals using this system
- Heights 50, 100, 200 mm.
- The drawers supported by a ISO drawer framework.
- Total opening, steel ball-bearing slides with retention system.
- Front stop, it protects the slides against strong impacts when closing.
- Card holder with identification cards ISO norm 7810.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>



Trolley prepared to mount on accessories at destination

Monobloc structure of AISI-304 stainless steel

Polyamide handle, PA, ergonomic  
Easy-to-clean

Polyethylene drawer front part, PE  
White, blue, lilac, green, yellow, red

**DOUBLE WHEELS**  
Ø 100 mm

2 Antistatic

2 with brake  
Ergonomic pedal

