

TECHNICAL FEATURES

100 SERIES SIMPLE TROLLEY

640 x 480 mm

EXTENSIVE MOBILY TROLLEYS

WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>



ISO DRAWERS 600 x 400 mm

- ABS ISO drawers of standard dimensions according to ISO 3394, compatible and interchangeable with ISO furniture of hospitals using this system
- Heights 50, 100, 200 mm.
- The drawers supported by a ISO drawer framework.
- Total opening, steel ball-bearing slides with retention system.
- Front stop, it protects the slides against strong impacts when closing.
- Card holder with identification cards ISO norm 7810.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>



Trolley prepared to mount on accessories at destination

Elevated push handle

Monobloc structure of AISI-304 stainless steel

Polyamide handle, PA, ergonomic
Easy-to-clean

Polyethylene drawer front part, PE
White, blue, lilac, green, yellow, red

SIMPLE WHEELS
Ø 125 mm

2 with brake

2 Antistatic



165 Y.AXWD12.WS12