

## TECHNICAL FEATURES

### 100 SERIES UNIT DOSE MEDICATION TROLLEY

640 x 480 mm

#### EXTENSIVE MOBILY TROLLEYS

#### WORK SURFACE

- Polyethylene, PE, high resistance: 100% pure alcohol (not diluted), bleach, detergent, soap solutions, oils, polish, aqueous ammonia, chlorhexidine.
- Protection on all 4 sides, including the front.

<http://www.insaustimedicaltrolleys.com/en/technical-features/work-surface/>





#### ISO DRAWERS 600 x 400 mm

- ABS of standard dimensions according to ISO 3394, interchangeable.
- Heights 50, 100, 200 mm.

<http://www.insaustimedicaltrolleys.com/en/technical-features/drawers/>

#### CASSETTE FOR MEDICATION, INTERCHANGEABLE, 5 PATIENTS

- Compatible with pharmacy ISO furniture.
- Polypropylene and aluminium grey colour.
- Cassette locked, epoxy painted steel, grey colour.
- Total weight 4.6 kg.



#### MEDICATION CONTAINER

- Capacity 3 litres.
- Identification card ISO norm 7810, international standard.
- Divider 1 unit included.

<https://www.insaustimedicaltrolleys.com/en/news/unit-dose-medication-trolleys/>

Trolley prepared to mount on accessories at destination

Elevated push handle

Polyamide handle, PA, ergonomic  
Easy-to-clean

Monobloc structure of AISI-304 stainless steel

Polyethylene drawer front part, PE  
White, blue, lilac, green, yellow, red

#### DOUBLE WHEELS

Ø 125 mm

2 Antistatic

2 with brake  
Ergonomic pedal

